



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

NO.330 Shichang Road, Weihai City, Shandong Province, 264209, China Phone +86 (0)631 5623761 | Mobile +86 183 6918 7836 | Fax +86 (0)631 5622563 Jameszhan@weigaogroup.com | exp_weigaogroup@163.com

DISPOSABLE BLOOD COLLECTION NEEDLE

KNBCSS Type Safety Blood Collection Needle


KNBCSS type safety blood collection needle is designed safety vein blood collection needle which connect safety scalp vein set with blood collection components. Push the safety protector after using up in order to shield off needle tip and avoid it stabbing nurse. The product have KNBCSS02 type and KNBCSS03 type with independent package for clinical choice.



OPERATION METHOD



1. Push the safety cover forward over the needle after using up



2. Make butterfly part totally lock into slot of safety device

KNBCSS02 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
25G	0.50	orange 🧶	25G×20	KNBCSS0225G20	3/4" (20mm)
23G	0.60	blue •	23Gx20	KNBCSS0223G20	3/4" (20mm)
22G	0.70	black •	22Gx20	KNBCSS0222G20	3/4" (20mm)
21G	0.80	green	21Gx20	KNBCSS0221G20	3/4" (20mm)

^{*}Remark: normal size in above table, we can also make special size as clients' request.

KNBCSS03 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
25G	0.50	orange 🥚	25Gx20	KNBCSS0325G20	3/4" (20mm)
23G	0.60	blue •	23Gx20	KNBCSS0323G20	3/4" (20mm)
22G	0.70	black •	22Gx20	KNBCSS0322G20	3/4" (20mm)
21G	0.80	green	21Gx20	KNBCSS0321G20	3/4" (20mm)


^{*}Remark: normal size in above table, we can also make special size as clients' request.

KNBCHS Type Safety Blood Collection Needle

KNBCHS type safety blood collection needle is designed safety vein blood collection needle which connect safety hypodermic needle with blood collection components. Push the safety protector after using up in order to shield off needle tip and avoid it stabbing nurse. This product have KNBCHS01、KNBCHS02 type and KNBCHS03 type with independent package for clinical choice.



OPERATION METHOD



- 1. After collection of the blood, holding thumb press down on the safety device of the holder.
- 2. Make safety device safety lock completely locked needle.



KNBCHS01 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
23G	0.60	Blue 🔵	23Gx25	KNBCHS0123G25	1" (25mm)
23G	0.60	Blue 🔵	23Gx32	KNBCHS0123G32	11/4" (32mm)
23G	0.60	Blue O	23Gx38	KNBCHS0123G38	11/2" (38mm)
22G	0.70	Black	22Gx38	KNBCHS0122G38	11/2" (38mm)
21G	0.80	Green 🔵	21Gx38	KNBCHS0121G38	11/2" (38mm)
20G	0.90	Yellow 👴	20Gx38	KNBCHS0120G38	11/2" (38mm)
18G	1.20	Pink 🌑	18Gx38	KNBCHS0118G38	11/2" (38mm)

^{*}Remark: normal size in above table, we can also make special size as clients' request.

KNBCHS02 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
23G	0.60	Blue O	23Gx25	KNBCHS0223G25	1" (25mm)
23G	0.60	Blue 🔵	23Gx32	KNBCHS0223G32	11/4" (32mm)
23G	0.60	Blue 🔵	23Gx38	KNBCHS0223G38	11/2" (38mm)
22G	0.70	Black	22Gx38	KNBCHS0222G38	11/2" (38mm)
21G	0.80	Green	21Gx38	KNBCHS0221G38	11/2" (38mm)
20G	0.90	Yellow 😑	20Gx38	KNBCHS0220G38	11/2" (38mm)
18G	1.20	Pink	18Gx38	KNBCHS0218G38	11/2" (38mm)

^{*}Remark: normal size in above table, we can also make special size as clients' request.

KNBCHS03 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
23G	0.60	Blue 🔵	23Gx25	KNBCHS0323G25	1" (25mm)
23G	0.60	Blue O	23Gx32	KNBCHS0323G32	11/4" (32mm)
23G	0.60	Blue O	23Gx38	KNBCHS0323G38	11/2" (38mm)
22G	0.70	Black	22Gx38	KNBCHS0322G38	11/2" (38mm)
21G	0.80	Green	21Gx38	KNBCHS0321G38	11/2" (38mm)
20G	0.90	Yellow 👴	20Gx38	KNBCHS0320G38	11/2" (38mm)
18G	1.20	Pink 🌑	18Gx38	KNBCHS0318G38	11/2" (38mm)

^{*}Remark: normal size in above table, we can also make special size as clients' request.

DISPOSABLE BLOOD COLLECTION NEEDLE

KNBCP Type Blood Collection Needle

Blood Collection Needle (Type KNBCP) is a comparatively economical needle. With its simple design and low inner volume of the tube, it can be easily handled and the blood sample can be fully collected by the clinic. Therefore, it is considered to be a very economical and useful blood collection needle.



KNBCP01 (Pen) TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
23G	0.60	blue 🔵	23Gx25	KNBCP0123G25	1" (25mm)
23G	0.60	blue 🔵	23Gx32	KNBCP0123G32	11/4" (32mm)
23G	0.60	blue 🔵	23Gx38	KNBCP0123G38	11/2" (38mm)
22G	0.70	black	22Gx38	KNBCP0122G38	11/2" (38mm)
21G	0.80	green 🌑	21Gx38	KNBCP0121G38	11/2" (38mm)
20G	0.90	yellow 😑	20Gx38	KNBCP0120G38	11/2" (38mm)
18G	1.20	pink 🌑	18Gx38	KNBCP0118G38	11/2" (38mm)

^{*}Remark: normal size in above table, we can also make special size as clients' request.

KNBCH01 (Injection Needle) TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
23G	0.60	blue 🔵	23Gx25	KNBCH 0123G25	1" (25mm)
23G	0.60	blue 🔵	23Gx32	KNBCH 0123G32	11/4" (32mm)
23G	0.60	blue 🔵	23Gx38	KNBCH 0123G38	11/2" (38mm)
22G	0.70	black	22Gx38	KNBCH 0122G38	11/2" (38mm)
21G	0.80	green 🌑	21Gx38	KNBCH 0121G38	11/2" (38mm)
20G	0.90	yellow 👴	20Gx38	KNBCH 0120G38	11/2" (38mm)
18G	1.20	pink 🔴	18Gx38	KNBCH 0118G38	11/2" (38mm)

^{*}Remark: normal size in above table, we can also make special size as clients' request.

KNBCS Type of Blood Collection Needle

Blood Collection Needle (Type KNBCS) is a medical instrument which combines the butterfly-winged infusion set with the blood collecting set. Now, it has been applied into clinical use. Through the end of the extension tube of blood collection needle, the transfusion condition can be well observed, which makes blood collection become successful.

Currently has produced single-winged blood collection needle (Type KNBCS01), butterfly-winged blood collection needle (Type KNBCS02) and butterfly-winged blood collection needle with holder (Type KNBCS03) for clinical uses.

KNBCS01 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
23G	0.60	blue 🔵	23Gx25	KNBCS0123G25	1" (25mm)
22G	0.70	black	22Gx25	KNBCS0122G25	1" (25mm)
21G	0.80	green 🔵	21Gx28	KNBCS0121G28	11/2" (38mm)

Remark: normal size in above table, we can also make special size as clients' request.

KNBCS02 TYPE

OD(mm)	Color	SIZE	REF	Length
0.50	orange 🧶	25Gx20	KNBCS0225G20	3/4" (20mm)
0.60	blue 🔵	23Gx20	KNBCS0223G20	3/4" (20mm)
0.70	black	22Gx20	KNBCS0222G20	3/4" (20mm)
0.80	green 🌑	21Gx20	KNBCS0221G20	3/4" (20mm)
0.90	yellow 🧶	20Gx20	KNBCS0220G20	3/4" (20mm)
	0.50 0.60 0.70 0.80	0.50 orange 0 0.60 blue 0 0.70 black 0 0.80 green 0	0.50 orange 25Gx20 0.60 blue 23Gx20 0.70 black 22Gx20 0.80 green 21Gx20	0.50 orange 25Gx20 KNBCS0225G20 0.60 blue 23Gx20 KNBCS0223G20 0.70 black 22Gx20 KNBCS0222G20 0.80 green 21Gx20 KNBCS0221G20

Remark: normal size in above table, we can also make special size as clients' request.

KNBCS03 TYPE

Gauge	OD(mm)	Color	SIZE	REF	Length
25G	0.50	orange 👴	25Gx20	KNBCS0325G20	3/4" (20mm)
23G	0.60	blue 🔵	23Gx20	KNBCS0323G20	3/4" (20mm)
22G	0.70	black	22Gx20	KNBCS0322G20	3/4" (20mm)
21G	0.80	green 🌑	21Gx20	KNBCS0321G20	3/4" (20mm)
20G	0.90	vellow 👴	20Gx20	KNBCS0320G20	3/4" (20mm)

Remark: normal size in above table, we can also make special size as clients' request.




KH Safety Holder

Blood Collection Needles with Safety Holder is a protective product applied in the clinical venous sampling which is achieved by blood collection needles to secure the operators. This product can be used in combination with the pen-type blood collection needle (Type KNBCP01) produced by our company and other blood collection needles of same structure in the market. And after the combination of two, the point of the needle can be perfectly isolated, which ensures that the needle will not accidentally hurt the doctors and patients. Thus, infection can be effectively reduced. Therefore, this product now has come to become the most common and safe auxiliary medical instrument of the blood collection needles in the application of the clinical department of hospital.



KH05 OPERATION METHOD



1. After collection of the blood, thumb press down on the safety device of the holder.




2. Make safety device safety lock completely locked needle.

KH06 OPERATION METHOD



1. Put blood collection needle into the holder.



Parts

BLOOD COLLECTION NEEDLE PARTS



NO.	Gauge	SIZE	REF	REMARK	DIAGRAM
1	20G	0.9x21mm	KNLA0220G21	AES	21 0
2	21G	0.8x19mm	KNLA0321G19	ABS	19 8
3	21G	0.8x19mm	KNLA0421G19	ABS	19
4	21G	0.8x19mm	KNLA0521G19	PP	19 8
5	21G	0.8x19mm	KNLA0621G19	ABS	
6	21G	0.8x19mm	KNLA0721G19	ABS	19 19

HOLDER PARTS FOR BLOOD COLLECTION NEEDLE



ARTWORK NO.	SIZE	REF	DIAGRAM
CTZ.41.00C	0.8×19	KNECH01	45 008 000 0
CTZ.176.00	0.9x19	KNECH02	47.8
CTZ.177.00	0.9x19	KNBCH03	47.8

HOLDER



ARTWORK NO.	REF	DIAGRAM
CTZ.lx.143	KH01	8 00 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0
CTZ.lx.153	KH02	9 E 99 0
CTZ.lx.152	KH06	47.8