



**BOTTLES
FRASCOS
FLACONS**

06



LABORATORY BOTTLE • FRASCO DE LABORATORIO • FLACON DE LABORATOIRE



*Blue screw cap, clear glass, autoclavable up to 140°C

*Tapón de rosca azul, vidrio blanco, esterilizable hasta 140°C

*Bouchon à vis bleu, verre transparent, autoclavable jusqu'à 140°C

BORO
3.3

ISO
4796

Retrace
Code

*Includes **pouring ring** for greater thermal and chemical resistance

*Incluye **anillo de vertido** para una mayor resistencia térmica y química.

*Comprend un **anneau verseur** pour une plus grande résistance thermique et chimique

Ref.	Cap. (ml)	GL	Pack Unit
21801145	25	25	10
21801175	50	32	10
21801245	100	45	10
21801365	250	45	10
21801445	500	45	10
21801545	1000	45	10
21801635	2000	45	10
21801735	5000	45	6
21801865	10000	45	1

LABORATORY BOTTLE • FRASCO DE LABORATORIO • FLACON DE LABORATOIRE



Red screw cap, clear glass, autoclavable up to 180°C

Tapón de rosca rojo, vidrio blanco, esterilizable hasta 180°C

Bouchon à vis rouge, verre transparent, autoclavable jusqu'à 180°C

BORO
3.3

ISO
4796

Retrace
Code

Ref.	Cap. (ml)	GL	Pack Unit
21802145	25	25	10
21802175	50	32	10
21802245	100	45	10
21802365	250	45	10
21802445	500	45	10
21802545	1000	45	10
21802635	2000	45	10
21802735	5000	45	6
21802865	10000	45	1

LABORATORY BOTTLE • FRASCO DE LABORATORIO • FLACON DE LABORATOIRE



Blue screw cap, amber glass, autoclavable up to 140°C

Tapón de rosca azul, vidrio ámbar, esterilizable hasta 140°C

Bouchon à vis bleu, verre ambre, autoclavable jusqu'à 140°C

BORO
3.3

ISO
4796

Retrace
Code

Ref.	Cap. (ml)	GL	Pack Unit
21804145	25	25	10
21804175	50	32	10
21804245	100	45	10
21804365	250	45	10
21804445	500	45	10
21804545	1000	45	10
21804635	2000	45	10
21804735	5000	45	6
21804865	10000	45	1



Retrace code
Código de rastreo
Code de traçabilité

BOD BOTTLE • FRASCO DBO • FLACON DBO

For incubating diluted samples of clean water, waste water, effluent water and polluted water during 5 day biochemical oxygen demand (BOD) test.

This test is used to determine the relative oxygen requirements of the tested sample.

Para la incubación de muestras de agua limpia, agua residual, agua efluente y contaminada durante el teste de 5 días de demanda biológica de oxígeno (DBO).

Este método pretende medir las necesidades relativas del oxígeno en la muestra en prueba.

Pour l'incubation des échantillons d'eau propre, eaux superficielles, eau usées et contaminées pendant l'essai des 5 jours pour la détermination de la demande biochimique en oxygène (DBO).

Cette méthode vise mesurer les besoins relatifs d'oxygène dans l'échantillon testé.

Ref.	Cap. (ml)	Pack Unit
9060067	300	1

→ Interchangeable stopper

BORO
3.3

→ Tapón intercambiable

→ Bouchon Interchangeable

**BOD BOTTLE • FRASCO DBO • FLACON DBO**

Ref.	IMG	Cap. (ml)	GL	Pack Unit
9060049	A	292	24/29	1
9060066	B	165	24/29	1

→ Interchangeable stopper

BORO
3.3

→ Tapón intercambiable

→ Bouchon Interchangeable

**WINKLER BOTTLE • FRASCO WINKLER • FLACON WINKLER**

Ref.	Cap. (ml)	Glass	Pack Unit
4471024	100	clear	1
4471036	250	clear	1
4471039	300	clear	1
4471136	250	amber	1

SODA

CALIBRATED



REAGENT BOTTLE • FRASCO REACTIVO • FLACON À RÉACTIFS



With stopper, narrow neck, clear glass

Con tapón, boca estrecha, vidrio blanco

Avec bouchon, col étroit, verre transparent

SODA

Ref.	Cap. (ml)	NS	Pack Unit
4116517	50	12/21	1
4116524	100	14/23	1
4116536	250	19/26	1
4116544	500	24/29	1
4116554	1000	29/32	1

REAGENT BOTTLE • FRASCO REACTIVO • FLACON À RÉACTIFS



With stopper, narrow neck, amber glass

Con tapón, boca estrecha, vidrio blanco

Avec bouchon, col étroit, verre transparent

SODA

Ref.	Cap. (ml)	NS	Pack Unit
4116817	50	12/21	1
4116824	100	14/23	1
4116836	250	19/26	1
4116844	500	24/29	1
4116854	1000	29/32	1

REAGENT BOTTLE • FRASCO REACTIVO • FLACON À RÉACTIFS



With stopper, wide neck, clear glass

Con tapón, boca ancha, vidrio blanco

Avec bouchon, col large, verre transparent

SODA

Ref.	Cap. (ml)	NS	Pack Unit
41185171	50	19/26	1
41185241	100	24/29	1
41185361	250	34/35	1
41185441	500	45/40	1
41185541	1000	45/40	1

REAGENT BOTTLE • FRASCO REACTIVO • FLACON À RÉACTIFS



With stopper, wide neck, amber glass

Con tapón, boca ancha, vidrio ámbar

Avec bouchon, col large, verre ambre

SODA

Ref.	Cap. (ml)	NS	Pack Unit
41188171	50	19/26	1
41188241	100	24/29	1
41188361	250	34/35	1
41188441	500	45/40	1
41188541	1000	45/40	1